





Appendix 2 - Compression options following a Lower Limb Assessment

Doncas

Clinical Commissioning Group

- Compression bandages are to be used in Tier 3/4 and District Nursing Services only.
- Compressionhosiery and wraps can be used in all Tiers and District Nursing Services following a Lower Limb Assessment being completed. • Please ensure the Lower Limb Assessment is in date as per Appendix 3.

Appropriate compression options depending on the clinical scenario						
Scenario	Hosiery Kit	Hosiery	Adjustable wraps	Compression bandages		
Normal leg shape	\checkmark	\checkmark	\checkmark	\checkmark		
Distortion due to oedema	×	×	\checkmark	🗸 Urgo K Two only		
Deep skin folds	×	Flat knit (ref to TVALS)	×	\checkmark		
Low to moderate exudate	\checkmark	\checkmark	\checkmark	\checkmark		
High exudate	×	×	×	\checkmark		
Self caring patient	\checkmark	\checkmark	\checkmark	×		
Carer involvement	\checkmark	\checkmark	\checkmark	×		

40mmHg (Full compression)	Hosiery Kit	Hosiery	Adjustable wraps	Compression bandages
	Altipress 40 or Activa Leg ulcer Kit (10mmHg liner with a 30mmHg over stocking)	Activa class 2 (18-24mmHg)	Juzo Wraps ACS Calf and Foot	Urgo K Two (18-25cm)
	Jobst Leg Ulcer Kit Liner or Medi Leg Ulcer Kit (20mmHg liner with a 20mmHg over stocking)	Actilymph class 2 (23-32mmHg)	Juxta Cures Calf (20mmHg - 40mmHg)	Urgo K Two (25-32cm)
		Sigvaris Essential Comfortable class 1 (18 - 21 mmHg) class 2 (23 - 32 mmHg)		Actico Wadding Actico Bandage 2 layer
		Mediven Plus class 2 (23-32mmHg)	Jobst Farrow wrap Lite (20-30 mmHg)	 4 layer system: Wadding Type 2 Support Bandage
		Jobst Bellavar/ Opaque/ Ambition for men/ Explore class 2 (23-32mmHg)	Jobst Farrow wrap strong (30-40mmHg)	Class 3 a BandageClass 3 b Bandage

	Hosiery Kit	Hosiery	Adjustable wraps	Compression bandages
20 mmHg (Reduced compression)	Medi Leg Ulcer Kit Liner (18- 21mmHg)	Activa class 1 (18-21mmHg)	Medi Juxta Lite (20mmHg - 40mmHg)	Urgo K Two Reduced (18-25cm)
	Jobst Leg Ulcer Kit Liner (18-21mmHg)	Actilymph class 1 (23-32mmHg)	Medi Juxta Cures (20mmHg - 40mmHg)	Urgo K Two Reduced (25-32cm)
		Sigvaris Opaque/Active class 1 (15-20mmHg)		Actico Actico Wadding Actico Bandage 1 layer
		Mediven Plus class 1 (18-21mmHg)		3 layer system:
		Jobst Opaque class 1 (18-21mmHg)	-	Wadding Type 2 Support Bandage
		Jobst Ambition for men/ Explore class 1 (18-21 mmHg)		TypeClass 3 a Bandage
	Hosiery Kit		Applicators	
10 mmHg (Liner)	Activa or Altipress Liner (10mmHg)	Actiglide		
		Rolly		

If the named product on this pathway is not available a temporary second line product is available to use. This can be found within the main text of the Doncaster Wide Wound Care Formulary Document.

Reference: Wound UK (2016) Best Practice Statement: Holistic Management of Venous Leg Ulceration.

Ezy- As (with or without handles)

Rolly

Developed by: The Tissue Viability and Lymphoedema Service and The Skin Integrity Team April 2021. For review June 2024.